

CoreTainer®

Stereotactic Core Biopsy Image and Transport System
With A/B Divider for Basket Collection

Protect Tissue Integrity: All-in-one system eliminates excessive handling of delicate core specimens by imaging and transporting in the same container.

Transport Safely: Accommodates formalin and is safe for transport to pathology with its watertight lid.

Image Clearly: Helps to eliminate potentially obscuring artifact caused by other makeshift methods. Micro-calcifications are clearly visible on the image.

Enhance Communication: Provides a section specifically for patient information, helping to reduce potential confusion by keeping cores and patient information together.

Save Time: Divided A/B compartments and radiolucent letters quickly allow correlation of core location to the image.

REF 303 **QTY** 12 / Box

- All-in-one system for separating, imaging, and transporting stereotactic core needle biopsy specimens with minimal handling
- Two separate compartments labeled A and B
- 3.5" diameter

Instructions for Use:


1. Use one CoreTainer® with A/B Divider for Basket Collection per biopsy site.
2. Complete patient information on CoreTainer label.
3. Place specimens into A and/or B compartments as desired.
4. Take image of specimens in CoreTainer*.
5. Indicate use of A and/or B compartments for pathology.
6. Lightly moisten specimens with saline or fix with formalin.
7. Seal lid on tightly before transporting CoreTainer to pathology.

*CoreTainer is designed for use with digital equipment. If foreign particles are introduced, they may appear during imaging.



Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com



Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

 Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010-7454 USA
Tel: 1.800.233.5539 or +1.860.583.4700
Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, CORETAINER, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and CORETAINER are Registered Community Trademarks. BEEKLEY is a registered trademark in Canada.

© 2023 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse 
- Non-Sterile 



If any serious incident occurs in relation to this device, report it to Beekley Medical and to the competent authority of the Member State in which the user and/or patient is established.

Beekley Medical®
when your diagnosis must be right™

CE **EC** **REP**
EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

IVD

REV: 303_INF-IFU-WEB_0423