

## When treatment ends, *the mark should too.*

**InkRight® semi-permanent tattoos** give radiation oncology teams a reliable way to mark patients for setup, without the lifelong reminder or discomfort of permanent tattoos. Care teams deliver treatment with confidence, and patients move forward without a lasting reminder of their care.

### Patient experience.

- No lasting reminder of treatment. No permanent ink.
- Pain-free application.
- Respects personal, cultural, and religious preferences.
- Suitable for patients of all ages.
- Stays visible throughout bathing and exercise.

### Clinical confidence.

- A reliable replacement for traditional tattoos and paint pens.
- Visible immediately and lasts up to two weeks.
- Easy realignment over existing marks for consistent, reproducible treatment setup for extended treatments.

### Safe and simple workflow.

- Needle-free application eliminates needle-stick risk.
- Individually packaged for single-patient use reducing risk of cross contamination.
- Peel and stick application.
- No water, scissors, or refrigeration needed.



**The InkRight Difference**

Patients increasingly object to permanent tattoos, especially in curative or short-course treatments where the marks linger long after care ends. InkRight removes that burden while preserving the clinical rigor radiation oncology demands.

**InkRight®**

SEMI-PERMANENT TATTOO  
| DOT 5.0



**REF 510**  
**QTY 50 Tattoos and 50 Primer Wipes / Box**  
**5.0mm diameter Semi-permanent Tattoo**

**InkRight®**

SEMI-PERMANENT TATTOO  
| CROSSHAIR



**REF 511**  
**QTY 50 Tattoos and 50 Primer Wipes / Box**  
**17mm x 17mm open crosshair Semi-permanent Tattoo**

THE OPEN CROSSHAIR ONLY  
FROM BEEKLEY MEDICAL

**InkRight®**

SEMI-PERMANENT TATTOO  
| DOT 2.5



**REF 512**  
**QTY 50 Tattoos and 50 Primer Wipes / Box**  
**2.5mm diameter Semi-permanent Tattoo**



**1.800.233.5539 • info@beekley.com • beekley.com**

Customers outside the U.S. – Contact your local distributor for pricing and product availability.  
To locate a distributor, call +1.860.583.4700 or email international@beekley.com

Manufactured by Beekley Corporation  
One Prestige Lane, Bristol, CT 06010-7454 USA  
Tel: 1.800.233.5539 or +1.860.583.4700  
Fax: 1.800.735.1234 beekley.com

Pat. – beekley.com/patents  
© – Reg. U.S. Pat. & Tm. Off.  
OPEN CROSSHAIR – Reg. U.S. Pat. & Tm. Off.  
© 2026 Beekley Corporation. All rights reserved.

**WARNINGS AND PRECAUTIONS:**

- Single Patient Use Only; Do Not Reuse ☒
- External Use Only
- Do Not Place on Compromised Skin
- Non-Sterile ⚠
- Do Not Use If Sensitive Or Allergic to Dyes Or Inks

If any serious incident occurs in relation to this device, report it to Beekley Medical and to the competent authority of the Member State in which the user and/or patient is established.



EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands



REV: IR-INF-0626