

PINPOINT®

Multi-Modality Fiducial Markers for CT / MRI Image Registration

fiducial markers with a unique conical design
to assist in accurate image registration

PinPoint provides...

ACCURACY


Unique conical design direct probes to 1.27mm diameter center hole—PinPoint's registration area is 37% smaller than other commonly used fiducial markers.

UNMISTAKABLE IMAGING

Bright and consistent imaging in both CT and MRI—PinPoint is easily identified on scans.

TIME SAVINGS

Easy application and removal—PinPoint's Pik-Up™ tab helps eliminate frustration and wasted time.



REF 128 QTY 30/Box

* product shown not actual size and not to scale.

Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010 USA
Tel: 1.800.233.5539 or +1.860.583.4700
Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, PINPOINT, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada. PIK-UP is an unregistered trademark of Beekley Corporation.

© 2019 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse (⊗)
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile (NON-STERILE)

CE EC REP EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



REV: PP128_INF_1019