

# TenderTouch™

## Cushioned Tape for Thermoplastic Masks

Beekley  
Medical®  
when your diagnosis  
must be right®

Protect patients from the rough edges of thermoplastic masks

Provide Your Patients with Maximum Comfort During Head and Neck Radiation Therapy Treatments with TenderTouch™ cushioned tape for thermoplastic masks

92% of Radiation Therapists stated that cutting the thermoplastic mask results in rough or sharp edges.<sup>1</sup>

### Comfort for the Patient; Ease of Use for the Radiation Therapist

# 1

#### Maximum Patient Comfort

Soft, cushioned tape covers ragged edges; relieves pressure points to provide comfort.

# 2

#### Stays on for Duration of Treatment

Sticks consistently to the thermoplastic material throughout the entire treatment period.

# 3

#### Durable Marking Surface

Fade-resistant; orientation marks remain visible.

# 4

#### Easy Application Saves Time

Flexible, scored-edged tape easily molds around edges; allowing for larger piece applications.

Perfect for masks that are cut because of claustrophobia, anxiety, tightness, bolus, bite blocks, electron fields, difficulty breathing, and/or sensitive skin.



#### Dosimetric Properties of TenderTouch™

The effects on surface dose of TenderTouch foam adhesive material is minimal and could be ignored in most clinical situations. In particular, the incremental change in surface dose for this material relative to the effect of the thermoplastic is not considered clinically significant in most foreseeable situations.

A copy of the full Dosimetrist Report is available upon request.

<sup>1</sup> Beekley Medical Market Survey July 18, 2018

# TenderTouch™

Cushioned Tape for Thermoplastic Masks




REF 418 QTY 200cm cut to measure roll / Box

2.5cm wide

**Place Your Order Today to Offer Patients Maximum Comfort**  
**1.800.233.5539 orders@beekley.com beekley.com**

Customers outside the U.S. – Contact your local distributor for pricing and product availability.  
To locate a distributor call +1.860.583.4700 or email [international@beekley.com](mailto:international@beekley.com)

Visit [beekley.com](http://beekley.com) for instructions for use.

 Manufactured by Beekley Corporation  
One Prestige Lane, Bristol, CT 06010-7454 USA  
Tel: 1.800.233.5539 or +1.860.583.4700 Fax: 1.800.735.1234 [beekley.com](http://beekley.com)

BEEKLEY, BEEKLEY MEDICAL, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada. TenderTouch is an unregistered trademark of Beekley Corporation. © 2023 Beekley Corporation. All rights reserved.

**WARNINGS AND PRECAUTIONS:**

- Single Patient Use Only; Do Not Reuse 
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Non-Sterile 
- Not made with natural rubber latex



EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands



REV: TT\_INF\_0523