


PINPOINT®

for small field of view imaging in MRI



PinPoint skin marker for small field of view imaging in MRI sits on top of the adhesive label and area of interest.



PinPoint skin marker for small field of view imaging in MRI provides a clean and unaltered view of the area of interest.

PRECISE

Images on 2 to 3 slices and assists technologists in communicating the EXACT area of interest.

EASY TO USE

Simply peel and stick. No more fumbling with tape or worrying about melting vitamin capsules.

UNMISTAKABLE

Bright and consistent imaging on every sequence allows for unmistakable communication between patient, technologist and radiologist.

ELIMINATES INDENTATION

PinPoint skin marker for small field of view imaging in MRI lies on top of a self-adhesive label, reducing indentation caused by old methods of marking.



6mm Radiance® filled pellet

REF 187 QTY 40/Box

Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010-7454 USA
Tel: 1.800.233.5539 or +1.860.583.4700
Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, PINPOINT, RADIANCE, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada. PIK-UP is an unregistered trademark of Beekley Corporation.

© 2023 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse (X)
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile (NON-STERILE)



EC REP

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



REV: PP187_INF_0723