


Best Practices with MR-SPOT® skin markers



Radiance® filled PinPoint® imaged on film


Choosing the Right Skin Marker for Small Field of View Imaging

When your patient says, “it hurts right here,” a skin marker designed to image on 2-3 slices will best communicate that exact area of interest to the interpreting radiologist.

See How PinPoint® skin markers for Small Field of View Compares to Larger Markers




Radiance® filled PinPoint skin marker for small field of view imaging on wrist




Radiance filled PinPoint skin marker for small field of view imaging imaged on film

Larger markers localize a much broader region than the true area of interest:



PinPoint® skin markers for small field of view imaging provides precise localization:



Beekley Medical® recommends PinPoint® skin marker [REF 187](#) for small field of view imaging

- Sized to be positioned precisely on a specific, identifiable area of interest
- Images on just 2 to 3 slices for accurate localization
- Bright consistent imaging on every sequence




PinPoint®

6mm Radiance® filled pellet

REF 187 QTY 40 / Box



Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

 Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010-7454 USA
Tel: 1.800.233.5539 or +1.860.583.4700 Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, MR-SPOT, PINPOINT, RADIANCE, SPOT, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and SPOT are Registered Community Trademarks. BEEKLEY is a registered trademark in Canada. PIK-UP is an unregistered trademark of Beekley Corporation.
© 2023 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse 
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile 



EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



REV: PP187_ISC_BP_WEB_0723