

All PointGuards® are hypoallergenic, latex-free devices that do not contain neither natural rubber latex nor dry natural rubber as a component in the product or its packaging.

The adhesive on PointGuards® is a hypoallergenic, medical-grade, skin contact pressure sensitive acrylic.

The use of the term "hypoallergenic" has come to indicate a product which is non-sensitizing to the general public. The hypoallergenic claim for this product is supported by clinical evaluation using the repeated insult patch test in humans, commonly known as the Draize test. To be termed hypoallergenic, products are required to show no evidence of sensitization potential under this test.

No subjects showed evidence of sensitization with PointGuards® adhesive in accordance with the ISO 10993 -1 "Biological Evaluation of Medical Devices", as put forth by the FDA.

Should you have any questions or need further assistance, please do not hesitate to contact us at 1.800.BEEKLEY (233.5539).