



Avoid misregistration of skin marks and the need for repeat needle adjustments with GuideLines® CT biopsy grids.


GuideLines®



Precision engineered and easy to use. Uniform spacing between parallel lines allows for single stick biopsies.



All-over medical grade adhesive contours to any body part without buckling – optimizing secure placement throughout the procedure.




REF 217 **QTY** 15 / Box

- 7" x 7.5" biopsy grid with 12 non-metallic radiopaque lines
- ideal for CT-guided biopsies of the abdomen and other large areas

THE GRID WITH RAISED LINES,
ONLY FROM BEEKLEY MEDICAL


Other Methods



Homemade grids are imprecise as well as time consuming to make and apply.



Other products tend to lift off patient making it difficult to generate an accurate mark in relation to the image.



REF 117 **QTY** 20 / Box

- 4" x 5" biopsy grid with 7 non-metallic radiopaque lines
- ideal for CT-guided biopsies of the neck and other small areas

THE GRID WITH RAISED LINES,
ONLY FROM BEEKLEY MEDICAL

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010-7454 USA
Tel: 1.800.233.5539 or +1.860.583.4700 Fax: 1.800.735.1234
beekley.com

BEEKLEY, BEEKLEY MEDICAL, GUIDELINES, THE GRID WITH RAISED LINES, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada. PIK-UP is an unregistered trademark of Beekley Corporation.

© 2023 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile



EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



REV: GL_INF_0523