


LOCALIZERS™

Diagnostic CT Skin Markers

Eliminate rescans caused by skin marker artifact.

Save time, reduce potential rescans, and increase patient throughput with Localizers – the first skin markers specifically calibrated for diagnostic CT

Localizers



- Image Brightly
- Low Density, Low Artifact
- Clean Image & Communication

Older Marking Methods



- Generate Considerable Artifact
- Can Obscure Areas of Concern
- May Cause Unnecessary Rescans



Localizers

REF 210 QTY 40 / Box

- 2.3mm low density pellet
- ideal for identifying soft tissue masses on the neck



Localizers

REF 211 QTY 40 / Box

- 4.0mm low density pellet
- ideal for identifying areas of interest on abdominal cases or with larger patients

Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010
Tel: 1.800.233.5539 or 1.860.583.4700
Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada. LOCALIZERS and PIK-UP are unregistered trademarks of Beekley Corporation.

© 2019 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



REV: LOC_INF_1019