

MR-SPOT[®]

**Beekley
Medical[®]**
when your diagnosis
must be right[®]

Medical professionals share how MR-SPOT[®] skin markers improved departmental communication and diagnostic accuracy in MRI

MRI Technologist, Hartford, CT

“Skin markers serve an important purpose: marking the anatomy so the doctors can see where there is a point of pain, a mass, or an area of interest. Without that marker, it can be hard to tell where something is hurting on an MRI. You really need to pinpoint it.

“I used to tape almonds, mayonnaise packs, or mustard packs to patients. They wouldn't show up well and when that happens, it's more work for the interpreting physician.

“I hated using the fish oil and vitamin E capsules. Imaging a total spine can take up to two hours. Those capsules would just dissolve and disappear. It would be nasty. Those capsules would melt and get all over the table and the patient. I knew there had to be something better.

“I have been using MR-SPOT skin markers for twenty years. They are so easy to use. Just place the markers and poof! They show up on every slice and the doctors are happy. They serve their purpose - marking the anatomy and pin-pointing so doctors can see the mass or point of pain.”

MRI Supervisor, New York, NY

“MR-SPOT skin markers help our technologists and radiologists communicate areas of focus, concern, and confirm coverage in MRI imaging. Vitamin capsules or fish oils do not appear bright on the image and can be easily missed. They can also melt over during exam and do not image after a few sequences.

“Our orthopedic hospital uses MR-SPOT skin markers for the long bone extremities to localize the area of pain. We also mark



areas like the ribs if there's pain after trauma or sports impact to ensure coverage on the MRI. Ortho-SPOT Packets[™] are used when we need a smaller skin marker for imaging hands, fingers, and feet or to denote a scar or tumor.”

MRI Supervisor, Evanston, IL

“We prefer MR-SPOT skin markers because they stick well. We use MR-SPOT Packets[™] to denote soft tissue masses, extremities, and spines as well as for breast and breast biopsy exams. These skin markers do not artifact in areas of interest. The markers are easy on the skin and patients are comfortable with them. Most importantly, they are bright on the image and the different sizes work well for all types of anatomy.”

Lead MRI Technologist, New Britain, CT

“It's frustrating when we rely on unprofessional skin markers like capsules; they don't serve their purpose. Doctors have a hard time discerning the area of interest and technologists that apply the markers are disappointed by the results. The vitamin E and fish oil capsules would not image clearly and you could only see them on a few sequences. They just weren't consistent or reliable.

“The (vitamin) capsules are pretty hard; not at all soft. Depending on what you're scanning and where that mass is that you're trying to localize, it could be very uncomfortable for the patient, especially if they have to lay there for a while. To



make matters more complicated, the capsules are kind of 'one-size fits all.' If I needed a longer marker I would have to make it myself and tape a couple of capsules together.

“MR-SPOT skin markers work very well! We apply them to the level of the thoracic spine we are scanning through, and we see them in all the images just like we need to; consistent and reliable.”

MR-SPOT®

**Beekley
Medical®**
when your diagnosis
must be right®

MR-SPOT® Skin Markers for MRI



PACKET MARKERS ONLY
FROM BEEKLEY MEDICAL

MR-SPOT Packets™

5.0cm Radiance® filled packet

REF 184 QTY 40 / Box

- for thoracic spines, larger areas of interest
- designed for optimum patient comfort



PACKET MARKERS ONLY
FROM BEEKLEY MEDICAL

MR-SPOT Packets™

1.75cm Radiance® filled packet

REF 185 QTY 40 / Box

- for extremities, soft tissue masses, and breast MRI
- designed for optimum patient comfort



PACKET MARKERS ONLY
FROM BEEKLEY MEDICAL

Ortho-SPOT Packets™

.75cm Radiance® filled packet

REF 186 QTY 40 / Box

- ideal for orthopedic, pediatric, or breast imaging in MRI
- designed for optimum patient comfort



MR-SPOT®

3.0cm Radiance® filled tube

REF 121 QTY 40 / Box

- for thoracic spines, larger areas of interest



MR-SPOT®

1.5cm Radiance® filled tube

REF 122 QTY 40 / Box

- for extremities, soft tissue masses, and breast MRI



PinPoint®

6mm Radiance® filled pellet

REF 187 QTY 40 / Box



- for small field of view imaging in MRI
- precisely identify an exact area of interest within MRI

**To learn more about MR-SPOT® Skin Markers for MRI, contact Beekley Medical®
at 1.800.233.5539 or email info@beekley.com**

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010-7454 USA
Tel: 1.800.233.5539 or +1.860.583.4700 Fax: 1.800.735.1234 beekley.com
BEEKLEY, BEEKLEY MEDICAL, MR-SPOT, ORTHO-SPOT, PINPOINT, RADIANCE, SPOT, the PACKET product configuration, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and SPOT are Registered Community Trademarks. BEEKLEY is a registered trademark in Canada. ORTHO-SPOT PACKETS, MR-SPOT PACKETS, and PIK-UP are unregistered trademarks of Beekley Corporation.
© 2023 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse 
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile 

   EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



REV: MRS_TM_0723