



# MR-SPOT®

Beekley's complete line of MRI skin markers —  
the ideal solution to all of your MRI skin marking needs


*Bright and consistent images on all sequences*



For smaller areas of concern



For minimal indentation



For sharpest, most distinct images



PACKET MARKERS ONLY  
FROM BEEKLEY MEDICAL

## Ortho-SPOT Packets™

.75cm Radiance® filled packet

**REF 186 QTY 40 / Box**

- for orthopedic imaging in MRI
- designed for precise localization and indentation-free imaging



PACKET MARKERS ONLY  
FROM BEEKLEY MEDICAL

## MR-SPOT Packets™

1.75cm Radiance® filled packet

**REF 185 QTY 40 / Box**

- for extremities, soft tissue masses
- designed for optimum patient comfort



PACKET MARKERS ONLY  
FROM BEEKLEY MEDICAL

## MR-SPOT Packets™

5.0cm Radiance® filled packet

**REF 184 QTY 40 / Box**

- for thoracic spines, larger areas of interest
- designed for optimum patient comfort



## PinPoint®

6mm Radiance® filled pellet

**REF 187 QTY 40 / Box**

- for small field of view imaging in MRI
- precisely identify an exact area of interest within MRI



## MR-SPOT®

1.5cm Radiance® filled tube

**REF 122 QTY 40 / Box**

- for extremities, soft tissue masses
- for sharp, distinct image on film



## MR-SPOT®

3.0cm Radiance® filled tube

**REF 121 QTY 40 / Box**

- for thoracic spines, larger areas of interest
- for sharp, distinct image on film

**Call 1.800.233.5539 • Fax 1.800.735.1234**  
**Visit beekley.com • Email info@beekley.com**

Customers outside the U.S. – Contact your local distributor for pricing and product availability.  
To locate a distributor call +1.860.583.4700 or email international@beekley.com



Manufactured by Beekley Corporation  
One Prestige Lane, Bristol, CT 06010 USA  
Tel: 1.800.233.5539 or +1.860.583.4700 Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, SPOT, PINPOINT, RADIANCE, MR-SPOT, ORTHO-SPOT, the PACKET product configuration, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and SPOT are Registered Community Trademarks. BEEKLEY is a registered trademark in Canada. MR-SPOT PACKETS, ORTHO-SPOT PACKETS, and PIK-UP are unregistered trademarks of Beekley Corporation.

© 2019 Beekley Corporation. All rights reserved.

#### WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile



EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



REV: MRS\_INF\_1019