



The Beekley Skin Marking System[®] for Mammography


A visual communication system of radiographically distinct skin markers.




Pellet = Nipple




Circle = Mole



Line = Surgical Scar



Triangle = Palpable Mass



Square = Area of Concern or Pain



beekley.com • 800.233.5539

Visit beekley.com for Product Safety Information

Beekley Medical[®]



Manufactured by Beekley Corporation, Bristol, CT 06010-7454 USA

BEEKLEY, BEEKLEY MEDICAL, BEEKLEY SKIN MARKING SYSTEM, the Raised Circle, the Raised Triangle, the Raised Square, and the Hummingbird design are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a registered Community Trademark. BEEKLEY is a registered trademark in Canada.
© 2023, 2025 Beekley Corporation. All rights reserved.