

ColdSPOT®

Disposable, Self-Adhesive Cold Therapy Packs

Comfortable, self-adhesive gel packs help reduce bleeding, bruising, pain, and inflammation following cyst aspirations, fine needle aspirations, minimally invasive breast biopsies.



ColdSPOT®

REF 500 QTY 100 / Case

- 2.75" x 3" gel pack
- self-adhesive side tabs
- for targeted and effective cold treatment following minimally invasive breast biopsy procedures



ColdSPOT® with Floral Design

REF 502 QTY 68 / Case

- 2" x 3" advanced gel pack
- self-adhesive side tabs
- for targeted and effective cold treatment following minimally invasive breast biopsy procedures



ColdSPOT® pillows

REF 505 QTY 50 / Box

- 2" x 3" advanced gel pack
- soft, flexible fabric covering
- self-adhesive side tabs
- for targeted and effective cold treatment following minimally invasive breast biopsy procedures

ColdSPOT® cold therapy packs have adhesive side tabs for exact positioning to eliminate shifting and the need for bulky bandages.

Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com



Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

Beekley Medical®
when your diagnosis must be right®

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010-7454 Tel: 1.800.233.5539 or +1.860.583.4700
Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, COLDSPOT, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada.
© 2023 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse 
- Recommended removal after no more than 2 hours to avoid potential adverse skin reaction to adhesive
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile 



EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



REV: CS_INF_0423