

Cradles®

Needle Localization Wire Protectors

Beekley Medical®
when your diagnosis
must be right®

Professional, Safe, Fast Wire Management

Cradles

- Clear, see-through
- Easy to apply
- Adheres and contours to skin
- Saves time over makeshift methods



Other Methods


- Hides wire
- Time consuming
- Unprofessional



Cradles

REF 309 QTY 12/Box

12" x 3" translucent cover with
foam adhesive border



Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010 USA
Tel: 1.800.233.5539 or +1.860.583.4700 Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, CRADLES and WHEN YOUR DIAGNOSIS MUST
BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark.
BEEKLEY is a registered trademark in Canada.

© 2019 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only;
Do Not Reuse (⊗)
- Remove After Use
- Non-Sterile (⚠)
- Not made with natural rubber latex



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



REV: CRA_INF_1119