

# Disposable X-Ray Markers for General Imaging

Communicate points of interest such as nipple locations, entry and exit sites, areas of pain, or anatomical orientation in general diagnostic and emergency room imaging.

## Right and Left X-Ray Markers:

### L-SPOT®

REF 112 QTY 64 / Box


- single-use, self-adhesive, disposable x-ray markers
- lead-free, radiopaque L for documenting left side on images
- sanitary single patient use



### R-SPOT®

REF 113 QTY 64 / Box

- single-use, self-adhesive, disposable x-ray markers
- lead-free, radiopaque R for documenting right side on images
- sanitary single patient use



## Skin Markers:

### D-SPOT®

REF 104 QTY 99 / Box

- lead-free radiopaque arrow with easy Pik-Up™ tab
- sanitary single patient use



### X-SPOT®

REF 101 QTY 150 / Box

- 1.5mm lead-free radiopaque pellet with easy Pik-Up™ tab
- sanitary single patient use



### Y-SPOT®

REF 102 QTY 132 / Box

- 2.0mm lead-free radiopaque pellet with easy Pik-Up™ tab
- sanitary single patient use



All markers are lead-free and not made with natural rubber latex.

Call 1.800.233.5539 • Fax 1.800.735.1234  
Visit [beekley.com](http://beekley.com) • Email [info@beekley.com](mailto:info@beekley.com)

Customers outside the U.S. – Contact your local distributor for pricing and product availability.  
To locate a distributor call +1.860.583.4700 or email [international@beekley.com](mailto:international@beekley.com)

Manufactured by Beekley Corporation

One Prestige Lane, Bristol, CT 06010  
Tel: 1.800.233.5539 or +1.860.583.4700  
Fax: 1.800.735.1234 [beekley.com](http://beekley.com)

BEEKLEY, BEEKLEY MEDICAL, SPOT, X-SPOT, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and SPOT are Registered Community Trademarks. BEEKLEY is a registered trademark in Canada. PIK-UP is an unregistered trademark of Beekley Corporation. © 2020 Beekley Corporation. All rights reserved.

#### WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse (X)
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile (NON STERILE)

CE EC REP EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

