

Disposable X-Ray Markers for General Imaging

Communicate points of interest such as nipple locations, entry and exit sites, areas of pain, or anatomical orientation in general diagnostic and emergency room imaging.

Right and Left X-Ray Markers:

L-SPOT®

REF 112 QTY 64 / Box


- single-use, self-adhesive, disposable x-ray markers
- lead-free, radiopaque L for documenting left side on images
- sanitary single patient use



R-SPOT®

REF 113 QTY 64 / Box

- single-use, self-adhesive, disposable x-ray markers
- lead-free, radiopaque R for documenting right side on images
- sanitary single patient use



Skin Markers:

D-SPOT®

REF 104 QTY 99 / Box

- lead-free radiopaque arrow with easy Pik-Up™ tab
- sanitary single patient use



X-SPOT®

REF 101 QTY 150 / Box

- 1.5mm lead-free radiopaque pellet with easy Pik-Up™ tab
- sanitary single patient use



Y-SPOT®

REF 102 QTY 132 / Box

- 2.0mm lead-free radiopaque pellet with easy Pik-Up™ tab
- sanitary single patient use



Call 1.800.233.5539 • Fax 1.800.735.1234
Visit beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1.860.583.4700 or email international@beekley.com

Manufactured by Beekley Corporation



One Prestige Lane, Bristol, CT 06010-7454

Tel: 1.800.233.5539 or +1.860.583.4700

Fax: 1.800.735.1234 beekley.com

BEEKLEY, BEEKLEY MEDICAL, L-SPOT, R-SPOT, SPOT, X-SPOT, Y-SPOT, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and SPOT are Registered Community Trademarks. BEEKLEY is a registered trademark in Canada. PIK-UP is an unregistered trademark of Beekley Corporation. © 2023 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse 
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile 

   EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



REV: GEN_INF_0523