




A NOVEL PATIENT STABILIZATION TECHNOLOGY



TECHNOLOGY



iFIX is a patient stabilization system comprised of both reusable and disposable components.

The disposable iFIX Fleece is attached to a reusable anchor product.



The unique micro anchors of the iFIX Patches and Adapters strongly secure the iFIX Fleece for efficient stabilization.



The comfortable stretchy fleece contours to the patient's anatomy like a second skin and secures the desired position.



Depending on your table's design, you can choose between the iFIX Patch that adheres to the table surface, the iFIX Slot Adapter that slides into the channels, or the iFIX Patch Adapter which connects to the fleece products.

IN JUST 4 SIMPLE STEPS THE PATIENT GETS STABILIZED



1a



1b

1

Estimate the required length of iFIX Fleece and cut.



2a



2b

2

On the opposite side of the table attach the iFIX Fleece to anchor product.

Anchor products:
iFIX Patch, iFIX Patch Adapter,
iFIX Slot Adapter



3a



3b

3

Pull iFIX Fleece tightly, up and over the patient.

Do not restrict patient's breathing.



4a



4b

4

Adhere the iFIX Fleece to the anchor product on the other side of the table and stabilize the patient.

If necessary, adjust iFIX Fleece to increase or decrease stabilization.

BENEFITS

simple
secure
intuitive

SAFE

head to toe
solution

VERSATILE

HYGENIC

disposable components lower risk for cross
contamination

EFFICIENT

no cleaning time
cost effective
fast user training

ifix

MULTI
MODALITY

supports several imaging
modalities
(CT, X-Ray, MRI, PET/SPECT, US)

Hip and lower spine stabilization



Hip and upper extremity stabilization



Head stabilization with the iFIX Head Mask



Head stabilization with head band



Lower extremity stabilization



Window for direct access



Distributed by:
Beekley Corporation
One Prestige Lane
Bristol, CT 06010 USA
Tel: 1-800-233-5539 or +1-860-583-4700

 **MANUFACTURED BY:**
iSYS Medizintechnik GmbH
Bergwerksweg 21
6370 Kitzbühel, Austria
Tel: +49 (8232) 99 41 956
info@interventional-systems.com

INTERVENTIONAL-SYSTEMS is a registered trademark
iSYS is a registered trademark
iFIX is a registered trademark
iFIX technology patent pending
© 2017 iSYS Medizintechnik GmbH

